

Product datasheet

Recombinant Human YY1 protein ab152809

1 Image

Description

Product name	Recombinant Human YY1 protein
Expression system	Wheat germ
Accession	P25490
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Human
Sequence	<p>MASGDTLYATDGSEMPAEIVELHEIEVETIPVETIETTIVGE EEEEDDD DEDGGGGDHGGGGGHGAGHHHHHHHHHHHPPMIALQP LVTDDPTQVHHH QEVILVQTREEVVGDDSDGLRAEDGFEDQILIPVPAPAG GDDDYIEQTL VTVAAGKSGGGGSSSSGGGRVKKGGGKKSGLKSYLSG GAGAAGGGGADP GNKKWEQKQVQIKTLEGEFSVTMWSSDEKKDIDHETVVE EQIIGENSPPD YSEYMTGKKLPPGGIPGIDLSDPKQLAEFARMKPRKIKED DAPRTIACPH KGCTKMFRDNSAMRKHLHHTHGPRVHVCAECGKAFVSS KLKRHQLVHTGE KPFQCTFEGCGKRFSLDFNLRTHVRIHTGDRPYVCPFDG CNKKFAQSTNL KSHILTHAKAKNNQ</p>
Predicted molecular weight	71 kDa including tags
Amino acids	1 to 414
Tags	GST tag N-Terminus

Specifications

Our [Abpromise guarantee](#) covers the use of **ab152809** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications ELISA

SDS-PAGE

Western blot

Form Liquid

Additional notes

Preparation and Storage

Stability and Storage Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
pH: 8.00
Constituents: 0.31% Glutathione, 0.79% Tris HCl


General Info

Function Multifunctional transcription factor that exhibits positive and negative control on a large number of cellular and viral genes by binding to sites overlapping the transcription start site. May play an important role in development and differentiation. The function of YY1 as an activator or a repressor is specified by the presence of other proteins. For example it acts as a repressor in absence of adenovirus E1A protein but as an activator in its presence.

Sequence similarities Belongs to the YY transcription factor family.
Contains 4 C2H2-type zinc fingers.

Cellular localization Nucleus matrix. Associated with the nuclear matrix.

Images



12.5% SDS-PAGE analysis of ab152809 stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors