

Product datasheet

Recombinant Human ZFP161 protein ab127611

Description

Product name	Recombinant Human ZFP161 protein
Expression system	Escherichia coli
Accession	O43829
Protein length	Protein fragment
Animal free	No
Nature	Recombinant
Species	Human
Predicted molecular weight	31 kDa
Amino acids	151 to 427
Tags	His-DHFR tag N-Terminus

Specifications

Our [Abpromise guarantee](#) covers the use of **ab127611** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE
Form	Lyophilized

Preparation and Storage

Stability and Storage	Shipped at 4°C. Store at -20°C. Constituents: 0.32% Tris HCl, 0.58% Sodium chloride
Reconstitution	Reconstitute with water to desired concentration.

General Info

Function	Transcriptional repressor of MYC and thymidine kinase promoters (By similarity). Transcriptional activator of the dopamine transporter (DAT), binding it's promoter at the consensus sequence 5'-CCTGCACAGTTCACGGA-3'.
Sequence similarities	Belongs to the krueppel C2H2-type zinc-finger protein family. Contains 1 BTB (POZ) domain.

	Contains 5 C2H2-type zinc fingers.
Domain	The BTB/POZ domain seems to direct the protein to discrete regions in the nucleus.
Cellular localization	Nucleus.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors