

Product datasheet

Recombinant Human Zinc finger protein 709 ab165753

1 Image

Overview

Product name	Recombinant Human Zinc finger protein 709
Protein length	Full length protein

Description

Nature	Recombinant
Source	Wheat germ
Amino Acid Sequence	
Accession	Q8N972
Species	Human

Sequence

MDSVVFEDVAVNFTQEEWALLGPSQKKLYRDVMQET
 FVNLASIGENWEEK
 NIEDHKNQGRKLRSHMVERLCERKEGSQFGETISQTP
 NPKPNKKTFRVK
 PYECVCGKDYMCHSSLNRHMRSHTEHRSYEHKYG
 EKS YECKE CKGRFS
 FRSSFRIHERTHTGEKPYKCKQCGKAFSWPSSFQIHE
 RTHTGEKPYECKE
 CGKAFYHTTFRGHMRMHTGEKPYKCKEKGKTFSHPS
 SFRNHERTHSGEK
 PYECKQCGKAFRYQTFQIHERTHTGEKPYQCKQCGK
 ALSCTSFRSHER
 IHTGEKPYKCKKCGKAFSFPSSFRKHERIHTGEKPYDC
 KECGKAFISLPS
 YRRHMIMHTGNGPYKCKEKGKAFDCPSSFQIHERTHT
 GEKPYECKQCGKA
 FSCSSSFRMHERTHTGEKPHECKQCGKAFSCSSSVR
 IHERTHTGEKPYEC
 KQCGKAFSCSSSFRMHERIHTGEKPYECKQCGKAFS
 FSSSFRMHERTHTG
 EKPYECKQCGKAFSCSSSFRMHERTHTGEKPYECKQ
 CGKAFSCSSSIRIH
 ERTHTGEKPYECKQCGKAFSCSSSVRMHERTHTGVK
 PYECKQCDKAFSCS
 RSFRIHERTHTGEKPYACQCGKAFKCSRSFRIHERVH

Amino acids	1 to 641
Tags	GST tag N-Terminus

Specifications

Our [Abpromise guarantee](#) covers the use of **ab165753** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	Western blot ELISA
Form	Liquid
Additional notes	Protein concentration is above or equal to 0.05 mg/ml.


Preparation and Storage

Stability and Storage	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.31% Glutathione, 0.79% Tris HCl
------------------------------	--

General Info

Function	May be involved in transcriptional regulation.
Sequence similarities	Belongs to the krueppel C2H2-type zinc-finger protein family. Contains 19 C2H2-type zinc fingers. Contains 1 KRAB domain.
Cellular localization	Nucleus.

Images



ab165753 on a 12.5% SDS-PAGE stained with Coomassie Blue.

SDS-PAGE - Recombinant Human Zinc finger protein 709 (ab165753)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors