

Product datasheet

Recombinant Human ZN528 protein ab164589

1 Image

Overview

Product name	Recombinant Human ZN528 protein
Protein length	Full length protein

Description

Nature	Recombinant
Source	Wheat germ
Amino Acid Sequence	
Species	Human

Sequence	<p>MALTQGPLKFMDVAIEFSQEEWKCLDPAQRTLYRDV MLENYRNLVSLGIC LPDLSVTSMLEQKRDPWTLQSEEKIANDPDGRECIKG VNTERRSSKLGNSA GNKPCKNQLGFTFQLHLSDLQLFQAERKISGCKHFEK PVSDNSSVSPLEK ISSSVKSHLLNKYRNNFDHAPLLPQEQKAHIREKAYKC NEHGQVFRASAS LTNQVIHNADNPYKCSECGKVFSCSSKLVHRRMHTGE KPYKCECGKLF SSNSNLSQHRIHTGEKPYKCECDKVFRRSSSKLAQH QRIHTGEKPYKCH ECDKVFNQIAHLVRHQKIHTGEKPYSCNKCGKVFSRH SYLAEHQTVHTGE KPYKCEECGKAFSVRSSLITHQLIHTGRKPYKCKEEDK VFGRKCFLTSHQ RIHTRERPYGCSQCCKIFSQKSDLIRHRKTHHTDEKPYK CNKCGTAFREFS DLTAHFLIHSGEKPYECKEKGKVFYKSSLTSHHRIHT GEKPYKCNRCGK VFSRSSNLVCHQKIHTGEKPYKCNQCGKVFNQASYLT RHQIHTGERPYR CSKCGKAFRGCSGLTAHLAITEKKSHECKECKGKIFTQ KSSLTNHHRIHI GEKPYKCTLCSKVFSHNSDLAQHQRVHS</p>
-----------------	--

Amino acids 1 to 628

Tags proprietary tag N-Terminus

Specifications

Our [Abpromise guarantee](#) covers the use of **ab164589** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	Western blot ELISA
Form	Liquid
Additional notes	Protein concentration is above or equal to 0.05 mg/ml.


Preparation and Storage

Stability and Storage	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.31% Glutathione, 0.79% Tris HCl
------------------------------	--

General Info

Function	May be involved in transcriptional regulation.
Sequence similarities	Belongs to the krueppel C2H2-type zinc-finger protein family. Contains 15 C2H2-type zinc fingers. Contains 1 KRAB domain.
Cellular localization	Nucleus.

Images



ab164589 on a 12.5% SDS-PAGE stained with Coomassie Blue.

SDS-PAGE - Recombinant Human ZN528 protein
(ab164589)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors