

## Product datasheet

# Recombinant Human ZNF133 protein ab127455

### Overview

---

<b>Product name</b>	Recombinant Human ZNF133 protein
<b>Protein length</b>	Protein fragment

### Description

---

<b>Nature</b>	Recombinant
<b>Source</b>	Escherichia coli
<b>Amino Acid Sequence</b>	
<b>Accession</b>	<a href="#">P52736-2</a>
<b>Species</b>	Human
<b>Molecular weight</b>	28 kDa
<b>Amino acids</b>	388 to 635
<b>Tags</b>	His-DHFR tag N-Terminus

### Specifications

---

Our [Abpromise guarantee](#) covers the use of **ab127455** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

<b>Applications</b>	SDS-PAGE
<b>Form</b>	Lyophilised

### Preparation and Storage

---

<b>Stability and Storage</b>	Shipped at 4°C. Store at -20°C. Constituents: 0.32% Tris HCl, 0.58% Sodium chloride
<b>Reconstitution</b>	Reconstitute with water to desired concentration.

### General Info

---

<b>Function</b>	May be involved in transcriptional regulation as a repressor.
<b>Tissue specificity</b>	Seems ubiquitous. Seen in the heart, brain, placenta, lung, liver, skeletal muscle, kidney and

pancreas.

**Sequence similarities**

Belongs to the krueppel C2H2-type zinc-finger protein family.  
Contains 15 C2H2-type zinc fingers.  
Contains 1 KRAB domain.

**Cellular localization**

Nucleus.

---

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

**Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

**Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors