

## Product datasheet

# Recombinant Human ZNF31 protein ab159839

1 Image

### Overview

<b>Product name</b>	Recombinant Human ZNF31 protein
<b>Protein length</b>	Full length protein

### Description

<b>Nature</b>	Recombinant
<b>Source</b>	Wheat germ
<b>Amino Acid Sequence</b>	
<b>Species</b>	Human
<b>Sequence</b>	<p>MAMALELQAQASPPQPEPEELLVKLEEDSWGSESKL            WEKDRGSVSGPEAS            RQRFRQFQYRDAAGPHEAFSQLWALCCRWLRPEIRL            KEQILELLVLEQFL            TILPREVQTTWVQARHPESGEEAVALVEDWHRETRTAG            QSGLELHTEETRP            LKTGEEAQSFQLQPVDPWPEGQSQKKGKNTCPDLP            NHLNAEVAPQPLKE            SGVPVSKPSNTSEKEQGPEFWGLSLINSKGRSTADYS            LDNEPAQALTWRD            SRAWEEQYQWDVEDMKVSGVHWGYEETKTFILAISE            SPFSEKLRTCHQNR            QVYRAIEQLRARGFLRTLQCRYRVKNLLRNRYKAKS            SHPPGTCPFYEE            LEALVRARTAIRATDGPGEAVALPRLGSDAEMDEQE            EGGWDPEEMAEDC            NGAGLVNVESTQGPRIAGAPALFQSRIGKNMGV</p>
<b>Amino acids</b>	1 to 433
<b>Tags</b>	proprietary tag N-Terminus

### Specifications

Our [Abpromise guarantee](#) covers the use of **ab159839** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

<b>Applications</b>	ELISA
---------------------	-------

	Western blot
<b>Form</b>	Liquid
<b>Additional notes</b>	Protein concentration is above or equal to 0.05 mg/ml.


## Preparation and Storage

<b>Stability and Storage</b>	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.31% Glutathione, 0.79% Tris HCl
------------------------------	--

## General Info

<b>Relevance</b>	ZSCAN20 belongs to the krueppel C2H2 type zinc finger protein family. It contains 10 C2H2 type zinc fingers and 1 SCAN box domain. ZSCAN20 may be involved in transcriptional regulation.
<b>Cellular localization</b>	Nucleus

## Images



ab159839 on a 12.5% SDS-PAGE stained with Coomassie Blue.

SDS-PAGE - Recombinant Human ZNF31 protein  
(ab159839)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors