

Product datasheet

Recombinant Human ZNF346 protein ab171689

1 Image

Overview

Product name	Recombinant Human ZNF346 protein
Protein length	Full length protein

Description

Nature	Recombinant
Source	Escherichia coli
Amino Acid Sequence	
Accession	Q9UL40
Species	Human

Sequence	<p>MGSSHHHHHH SSGLVPRGSH MGSMEYPAPA TVQAADGGAA GPYSSSELLE GQEPDGVRFD RERARRLWEA VSGAQPVGRE EVEHMIQKNQ CLFTNTQCKV CCALLISESQ KLAHYQSKKH ANKVKRYLAI HGMETLKGET KKLDSQKSS RSKDKNQCCP ICNMTFSSPV VAQSHYLGKT HAKNLKQKQ STKVEALHQN REMIDPDKFC SLCHATFNDP VMAQQHYVGK KHRKQETKLG LMARYGRLAD PAVTDFPAGK GYPCKTCKIV LNSIEQYQAH VSGFKHKNQS PKTVASSLGQ IPMQRQPIQK DSTTLED</p>
-----------------	---

Molecular weight	35 kDa including tags
Amino acids	1 to 294
Tags	His tag N-Terminus

Specifications

Our [Abpromise guarantee](#) covers the use of **ab171689** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE
	Mass Spectrometry
Mass spectrometry	MALDI-TOF

Purity >90% by SDS-PAGE.
ab171689 is purified using conventional chromatography techniques.

Form Liquid

Preparation and Storage


Stability and Storage Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituents: 0.03% DTT, 0.32% Tris HCl, 50% Glycerol, 1.17% Sodium chloride

General Info

Images



15% SDS-PAGE analysis of ab171689 (3µg).

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors