

Product datasheet

Recombinant Human ZNF408 protein ab164114

1 Image

Overview

Product name	Recombinant Human ZNF408 protein
Protein length	Full length protein

Description

Nature	Recombinant
Source	Wheat germ
Amino Acid Sequence	
Species	Human

Sequence	<p>MEEAEELLLEGKKALQLAREPRLGLDLGWNPSGEGC TQGLKDVPPEPTRD ILALKSLPRGLALGPSLAKEQRLGVWCVDPLQPGLL WGPLEEEESASKEK GEGVKPRQEENLSLGPWGDVCACEQSSGWTSLVQR GRLESEGNVAPVRIS ERLHLQVYQLVLPGSELLLPQPSSEGPSLTQPGLDK EAAVAVVTEVESA VQQEVASPGEDAAEPCIDPGSQSPSGIQAENMVSPG LKFPDQDRISKDSQ PLGPLLQDGDVDEECPAQAQMPPPELQNSATQQDP DGSGASFSSSARGTQ PHGYLAKKLHSPSDQCP PRAKTPEPGAQQSGFPTLS RSPPGPAGSSPKQG RRYRCGECGKAFLQLCHLKKHAFVHTGHKPFLCTECG KSYSSEESFKAHM LGHRGVRPFPCPCDKAYGTQRDLKEHQVVHSGARP FACDQCGKAFARRP SLRLHRKTHQVPAAPAPCPCPVCGRPLANQGSLRNH MRLHTGEKPFLCPH CGRAFRQRGNLRGHLRLHTGERPYRCPHCADAFPQL PELRRHLISHTGEA HLCPVCGKALRDPHTLRAHERLHSGERPFPCPCQCGR AYTLATKLRRHLKS HLEDKPYRCPTCGMGYTLPQSLRRHQLSHRPEAPCS PPSVPSAASEPTVV</p>
-----------------	---

LLQAEPQLLDTHREEEVSPARDVVEVTISESQEKCFV
VPEEPDAAPSLVL IHKDMGLGAWAEVVEVEMGT

Amino acids	1 to 720
Tags	proprietary tag N-Terminus

Specifications

Our [Abpromise guarantee](#) covers the use of **ab164114** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	Western blot ELISA
Form	Liquid
Additional notes	Protein concentration is above or equal to 0.05 mg/ml.


Preparation and Storage

Stability and Storage	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.31% Glutathione, 0.79% Tris HCl
------------------------------	--

General Info

Function	May be involved in transcriptional regulation.
Sequence similarities	Contains 10 C2H2-type zinc fingers.
Cellular localization	Nucleus.

Images



ab164114 on a 12.5% SDS-PAGE stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors