

Product datasheet

Recombinant Human ZNF658 protein ab161909

1 Image

Overview

Product name	Recombinant Human ZNF658 protein
Protein length	Full length protein

Description

Nature	Recombinant
Source	Wheat germ
Amino Acid Sequence	
Species	Human

Sequence	<p>MSQASVSFQDVTVEFTREEWQH LGPVERTLYRDV ML ENYSHLISVGYCIT KPKVSKLEKGEEPWSLEDEF LNQRYPGYFKVDHIKGI REKQEKLWQEI FISDADKTL SKEGQKVLEKPFNLEIAPELSEKISCKCDS HRMNLPVASQL IISERKYSRKKTEYMNVCEKLQLDIKHEKAHAEKSYEH GENAKAFSYKK DQHWKFQTL EESFECDGSGQGLYDKTICITPQSFLTGE KSCKDDEF RKNF DKITL FNHMRTDTRGKCS DLNEYGTSCDKTTAVEYNKV HMAMTHYECNER GINFSRKSPLTQSQR TITGWSAFESNKCEENFSQSSA HIVHQKTQAGDKF GEHNECTDALYQKLDFTA HQRIHTEDKFYLSDEHGKC RKSFYRKAHLIQH QRPHSGEKTYQYEECAK SFCSSSHPIQHPGTYVGFKL YECNECGKAFCQN SNLSKHLRIHTKEKPCD NNGCGRSYKSPLIGHQKTDAE MELCGGSEYGKT SHLKGHQ RILMGEKPYECIECGKTF SKTSHLRAHQRIHT GEKPYECVECE KTFSHKTHLSVHQRVHTGEKPYECNDCGKSFTYNSAL RAHQRIHTGLAVC HVGILMATKSLFHQS</p>
-----------------	--

Amino acids	1 to 615
--------------------	----------

Tags proprietary tag N-Terminus

Specifications

Our [Abpromise guarantee](#) covers the use of **ab161909** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications ELISA

Western blot

Form Liquid

Additional notes Protein concentration is above or equal to 0.05 mg/ml.

Preparation and Storage

Stability and Storage Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00

Constituents: 0.31% Glutathione, 0.79% Tris HCl

General Info

Function May be involved in transcriptional regulation.


Sequence similarities Belongs to the krueppel C2H2-type zinc-finger protein family.

Contains 24 C2H2-type zinc fingers.

Contains 1 KRAB domain.

Cellular localization Nucleus.

Images



ab161909 on a 12.5% SDS-PAGE stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors