

Product datasheet

Recombinant Human ZNF688 protein ab165537

1 Image

Overview

Product name	Recombinant Human ZNF688 protein
Protein length	Full length protein

Description

Nature	Recombinant
Source	Wheat germ
Amino Acid Sequence	
Species	Human
Sequence	<p>MAPPPAPLLAPRPGETRPGCRKPGTVSFADVAVYFS PEEWGCLRPAQRAL YRDVMQETYGHLGALGFPGPKPALISWMEQESEAWS PAAQDPEKGERLGG ARRGDVPNRKEEEPEEVPRAKGPRKAPVKESPEVLV ERNPDPAISVAPAR AQPPKNAAWDPTTGAQPPAPIPSMDAQAGQRRHVCT DCGRRFTYPSLLVS HRRMHSGERPFPCPECGMRFKRKF AVEAHQWIHRSC SGGRRGRRPGIRAV PRAPVRGDRDPPVLFRRHYPDIFECCG</p>
Amino acids	1 to 276
Tags	GST tag N-Terminus

Specifications

Our [Abpromise guarantee](#) covers the use of **ab165537** in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	ELISA Western blot
Form	Liquid
Additional notes	Protein concentration is above or equal to 0.05 mg/ml.

Preparation and Storage

Stability and Storage

Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00

Constituents: 0.31% Glutathione, 0.79% Tris HCl

General Info

Function

May be involved in transcriptional regulation.

Sequence similarities

Belongs to the krueppel C2H2-type zinc-finger protein family.


Contains 2 C2H2-type zinc fingers.

Contains 1 KRAB domain.

Cellular localization

Nucleus.

Images



SDS-PAGE - Recombinant Human ZNF688 protein
(ab165537)

ab165537 on a 12.5% SDS-PAGE stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors