

Product datasheet

Recombinant Influenza A Nucleoprotein (His tag)
ab214834

1 Image

Description

Product name	Recombinant Influenza A Nucleoprotein (His tag)	
Purity	> 95 % SDS-PAGE.	
Expression system	Escherichia coli	
Accession	B6A6U5	
Protein length	Full length protein	
Animal free	No	
Nature	Recombinant	
Species	Influenza A	
Sequence	<pre> MASQGTKRSYEQMETDGDQRNATEIRASVGKMIDGIG RFYIQMCTELKLS DHEGRLIQNSLTIEKMVLSAFDERRNKYLEEHPSAGKD PKKTGGPIYRRV DGKWMRELVLVDKKEIRRW RQANNGEDATSGLTHIMWHSNLNDATY QRTRALVRTGMDPRMCSLMQGSTLPRRSGAAGAAVK GIGTMVMELIRMVK RGINDRNFWRGENGRKTRSAYERMCNILKGKFQTAQ RAMVD QVRESR NPGNAEIEDLIFLARSALILRGSVAHKSCLPACAYGPAV SSGYDFEKEGY SLVGIDPFKLLQNSQIYSLIRPNENPAHKSQLVWMACH SAAFEDLRLLSF IRGTKVSPRGKLST RGVQIASNENMDNMGSSTLELRSGYWAI RTRSGG NTNQQRASAGQTSVQPTFSVQRNLPF EKSTIMAAFTGNTEGRTSDMRA EIIRMMEGAKPEEVSFRRGRGVFELSDEKATNPVPSFD MSNEGSYFFGDN AEEYDN </pre>	
Predicted molecular weight	57 kDa including tags	
Amino acids	1 to 498	
Tags	His tag C-Terminus	
Additional sequence information	(A/Brisbane/10/2007) (H3N2) (Genebank No. CY035025).	

Specifications

Our [Abpromise guarantee](#) covers the use of **ab214834** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications SDS-PAGE


Form Liquid

Preparation and Storage

Stability and Storage Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.

Constituents: PBS, 0.1% BSA, 25% Glycerol

Images



4-20% SDS-PAGE analysis of ab214834.

SDS-PAGE - Recombinant Influenza A
Nucleoprotein (His tag) (ab214834)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors