

Product datasheet

Recombinant Influenza A virus (H5N1/HA1) protein (His tag) ab217654

1 Image

Description

<b>Product name</b>	Recombinant Influenza A virus (H5N1/HA1) protein (His tag)
<b>Purity</b>	> 95 % SDS-PAGE.
<b>Expression system</b>	HEK 293 cells
<b>Accession</b>	<a href="#">C4P278</a>
<b>Protein length</b>	Protein fragment
<b>Animal free</b>	No
<b>Nature</b>	Recombinant
<b>Sequence</b>	<p>QIC IGYHANNSTE QVDTIMEKNV TVTHAQDILE            KTHNGKLCNL DGVKPLILKD CSVAGWLLGN            PMCDEFLNVS EWSYVEKAS PANGLCYPGD            FNDYEELKHL LSRINHLKKI KIIPKSYWSN HEASSGVSA            CSYLGEPSFF RNVVWLIKKN NTYPPIKVNY TNTNQEDLLV            LWGIHHPNDE KEQIRYQNP NTSISVGTST LNQLVPKIA            TRPKVNGQSG RMEFFWTILK PNDSINFDSN            GNFIAPYAY KIAKKGDSVI MKSELEYGNC            NTKCQTPMGA INSSMPFHNI HPLTIGECPK YVKS NRLVLA            TGLRNAPQRE GGRRRKRGLF GAIAGFIEGG            WQGMVDGWYG YHHSNEQGSG YAADKESTQK            AIDGITNKVN SIIDKMNTQF EIVGREFNNL ERRIENLNKK            MEDGFLDVWT YNAELLV LME NERTLDFHDS            NVKNLYEKVR LQLRDN AKEL GNGCFEFYHK            CDNECMESVR NGTYDYPQYS EEARLSREEI            SGVKMESMVT YQIL</p>
<b>Amino acids</b>	18 to 534
<b>Tags</b>	His tag C-Terminus
<b>Additional sequence information</b>	Influenza A (A/chicken/Vietnam/NCVD-016/2008). (Genebank No. ACO07033).

Specifications

Our [Abpromise guarantee](#) covers the use of **ab217654** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

**Applications** SDS-PAGE

**Form** Liquid

### Preparation and Storage


---

**Stability and Storage** Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C long term.

Constituent: 100% PBS

### Images

---



SDS-PAGE analysis of ab217654.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors