

Product datasheet

Recombinant Intimin protein (His tag) ab225566

1 Image

Description

| | |
|--|---|
| Product name | Recombinant Intimin protein (His tag) |
| Purity | > 90 % SDS-PAGE. |
| Expression system | Escherichia coli |
| Accession | P52869 |
| Protein length | Protein fragment |
| Animal free | No |
| Nature | Recombinant |
| Species | Hafnia alvei |
| Sequence | <pre> ITEIKADKTTAVANGKDAVTTYTKVMKDGKPLSGEEVFTTT TLGTLKST EKTNTNGYRKVSLTSANQGKSLVSASVTMPQLMLKLELV EFFTQLTIDNG NVEIVGTGAKGKLPNVWLQYGQVNLKANGGNGKYTWYSA NPAIASVDPSS GQVTLKDKGETTITVVS GDKQTAIYTIAMPNSIVSVNSSGRV DYNTANNI CKNIKGSLPSSIKELKDLYDDWGAANKYQHYSQESITAWTL QTSENKVQG VASTYDLVRKNPLIDKVDIAGNYAYAVCVK </pre> |
| Predicted molecular weight | 46 kDa including tags |
| Amino acids | 1 to 280 |
| Tags | His tag N-Terminus |
| Additional sequence information | 6xHis-SUMO tag at the N-terminus. |

Specifications

Our [Abpromise guarantee](#) covers the use of **ab225566** in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| | |
|---------------------|----------|
| Applications | SDS-PAGE |
| Form | Liquid |

Preparation and Storage

Stability and Storage


Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
Constituents: 50% Glycerol (glycerin, glycerine), Tris buffer

General Info

Cellular localization

Cell outer membrane

Images



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) analysis of ab225566 with 5% enrichment gel and 15% separation gel.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors