

Product datasheet

Recombinant Measles virus Priorix, Schwarz strain, nucleocapsid protein ab74559

3 References 1 Image

Description

<b>Product name</b>	Recombinant Measles virus Priorix, Schwarz strain, nucleocapsid protein	
<b>Purity</b>	> 95 % SDS-PAGE. ab74559 is purified by ultracentrifugation.	
<b>Expression system</b>	Saccharomyces cerevisiae	
<b>Accession</b>	<a href="#">AAF85699.1</a>	
<b>Protein length</b>	Full length protein	
<b>Animal free</b>	No	
<b>Nature</b>	Recombinant	
<b>Sequence</b>	<p>MATLLRSLALFKRNKDKPPITSGSGGAIRGIKHIIVPIPGDSS ITTRSR LLDRLVRLIGNPDVSGPKLTGALIGLSLFFVESPGQLIQRITD DPDVSIR LLEVVSQSDQSQSGLTFASRGTNMEDEADQYFSDDDPISS DQSRFGWFGNK EISDIEVQDPEGFNMILGTILAQIWWLLAKAVTAPDTAADSE LRRWIKYT QQRVVGEFRLERKWL DVVRNRIAEDLSLRRFMVALILDIK RTPGNKPRI AEMICDIDTYVEAGLASFILTKFGIETMYPALGLHEFAGELS TLESLM NLYQQMGETAPYMVILENSIQNKFSAGSYPLLWSYAMGVG VELENSMGGL NFGRSYFDPAYFRLGQEMVRRSAGKVSSTLASELGITAED ARLVSEIAMH TTEDKISRAGVPRQAQVSFLHGDQSENELPRLGGKEDRR VKQSRGEARES YRETGPSRASDARAAHLPTGTPLDIDTATESSQDPQDSRR SADALLRLQA MAGISEEQGSDTDTPIVYNDNRLLD</p>	
<b>Predicted molecular weight</b>	58 kDa	
<b>Amino acids</b>	1 to 525	

Specifications

Our [Abpromise guarantee](#) covers the use of **ab74559** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

<b>Applications</b>	Western blot ELISA SDS-PAGE
<b>Form</b>	Lyophilized

---

## Preparation and Storage

<b>Stability and Storage</b>	Shipped at 4°C. After reconstitution store at -20°C. Avoid freeze / thaw cycles. Constituent: PBS
<b>Reconstitution</b>	Lyophilized, reconstitute with 120µl deionised water and add glycerol to 50%


---

## General Info

<b>Relevance</b>	Priorix is a new investigational measles, mumps and rubella (MMR) vaccine containing the Schwarz measles strain, the RA 27/3 rubella strain and a new mumps strain RIT 4385.
------------------	--

---

## Images



SDS-PAGE showing ab74559 at approximately 60kDa. (2 µg/lane).

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

---

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors