

Product datasheet

Recombinant mouse PDGF BB protein ab78593

2 Images

Description

Product name	Recombinant mouse PDGF BB protein
Purity	= 11 % SDS-PAGE. Purity Greater than 98% as determined by: Reducing and non-reducing SDS-PAGE. Analytical HPLC. Endotoxin level as measured by LAL is <0.01ng/ug or <0.1EU/ug. This product was sterile filtered.
Expression system	Escherichia coli
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Mouse
Sequence	MSLGSLAAAE PAVIAECKTR TEVFQISRNL IDRTNANFLV WPPCVEVQRC SGCCNNRNVQ CRASQVQMRP VQVRKIEVR KKPIFKKATV TLEDHLACKC ETIVTPRPVT

Specifications

Our [Abpromise guarantee](#) covers the use of **ab78593** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE
Form	Lyophilized
Additional notes	It is recommended to add a carrier protein (0.1% HSA or BSA) for long term storage.


Preparation and Storage

Stability and Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. Constituent: 0.294% Sodium citrate This product is an active protein and may elicit a biological response in vivo, handle with caution.
Reconstitution	Reconstitute the lyophilized product with sterile water at a concentration of 0.1 - 0.5 mg's/ml, which can be further diluted into other aqueous solutions.

Images



SDS PAGE analysis of ab78593 under non-reducing (-) and reducing (+) conditions. Stained with Coomassie Blue.



ab78593 used in Functional Studies.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors