

Product datasheet

Recombinant Mouse PPAR gamma protein (His tag)  
ab236330

1 Image

Description

<b>Product name</b>	Recombinant Mouse PPAR gamma protein (His tag)
<b>Purity</b>	> 85 % SDS-PAGE.
<b>Expression system</b>	Escherichia coli
<b>Accession</b>	<a href="#">P37238</a>
<b>Protein length</b>	Full length protein
<b>Animal free</b>	No
<b>Nature</b>	Recombinant
<b>Species</b>	Mouse
<b>Sequence</b>	<p>MGETLGDSPVDPEHGAFADALPMSTSQEITMVDTEMPF  WPTNFGISSVDL  SVMEDHSHSFDIKPFTTVDFSSISAPHYEDIPFTRADPMVA  DYKYDLKLQ  EYQSAIKVEPASPPYYSEKTQLYNRPHEEPSNSLMAIECR  VCGDKASGFH  YGVHACEGCKGFFRRTIRLKLIDRCDLNCRHKSRNKCC  YCRFQKCLA  VGMSHNAIRFGRMPQAEKEKLLAEISSDIDQLNPESADLR  ALAKHLYDSY  IKSFPLTKAKARAILTGKTTDKSPFVIYDMNSLMMGEDKIKF  KHITPLQE  QSKEVAIRIFQGCQFRSVEAVQEITEYAKNIPGFINLDLNDQ  VTLLKYGV  HEIYTMLASLMNKDGVLISEGQGFMTREFLKSLRKPFGDF  MEPKFEFAV  KFNALELDDSDLAIFIAVILSGDRPGLLNKPIEDIQDNLLQ  ALELQLK  LNHPESQLFAKVLQKMTDLRQVTEHVQLLHVIKKTETDM  SLHPLLQEIYKDLY</p>
<b>Predicted molecular weight</b>	63 kDa including tags
<b>Amino acids</b>	1 to 505
<b>Tags</b>	His tag N-Terminus

**Additional sequence information** N-terminal 10xHis-tagged.

## Specifications

---

Our [Abpromise guarantee](#) covers the use of **ab236330** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

**Applications** SDS-PAGE

**Form** Liquid

## Preparation and Storage

---

**Stability and Storage** Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.  
Constituents: Tris buffer, 50% Glycerol (glycerin, glycerine)

## General Info

---

**Function** Receptor that binds peroxisome proliferators such as hypolipidemic drugs and fatty acids. Once activated by a ligand, the receptor binds to a promoter element in the gene for acyl-CoA oxidase and activates its transcription. It therefore controls the peroxisomal beta-oxidation pathway of fatty acids. Key regulator of adipocyte differentiation and glucose homeostasis.

**Tissue specificity** Highest expression in adipose tissue. Lower in skeletal muscle, spleen, heart and liver. Also detectable in placenta, lung and ovary.

**Involvement in disease** Note=Defects in PPARG can lead to type 2 insulin-resistant diabetes and hypertension. PPARG mutations may be associated with colon cancer.  
Defects in PPARG may be associated with susceptibility to obesity (OBESITY) [MIM:601665]. It is a condition characterized by an increase of body weight beyond the limitation of skeletal and physical requirements, as the result of excessive accumulation of body fat.  
Defects in PPARG are the cause of familial partial lipodystrophy type 3 (FPLD3) [MIM:604367]. Familial partial lipodystrophies (FPLD) are a heterogeneous group of genetic disorders characterized by marked loss of subcutaneous (sc) fat from the extremities. Affected individuals show an increased preponderance of insulin resistance, diabetes mellitus and dyslipidemia.  
Genetic variations in PPARG can be associated with susceptibility to glioma type 1 (GLM1) [MIM:137800]. Gliomas are central nervous system neoplasms derived from glial cells and comprise astrocytomas, glioblastoma multiforme, oligodendrogliomas, and ependymomas.  
Note=Polymorphic PPARG alleles have been found to be significantly over-represented among a cohort of American patients with sporadic glioblastoma multiforme suggesting a possible contribution to disease susceptibility.


**Sequence similarities** Belongs to the nuclear hormone receptor family. NR1 subfamily.  
Contains 1 nuclear receptor DNA-binding domain.

**Cellular localization** Nucleus.

---

## Images

---



SDS-PAGE - Recombinant Mouse PPAR gamma protein (His tag) (ab236330)

(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) analysis with 5% enrichment gel and 15% separation gel of ab236330.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors