

Product datasheet

Recombinant Mouse SERPINE2/PN-1 protein (His tag)
ab222348

1 Image

Description

Product name	Recombinant Mouse SERPINE2/PN-1 protein (His tag)	
Purity	> 95 % SDS-PAGE. ab222348 was purified using conventional chromatography techniques.	
Endotoxin level	< 1.000 Eu/µg	
Expression system	Baculovirus infected insect cells	
Accession	Q07235	
Protein length	Full length protein	
Animal free	No	
Nature	Recombinant	
Species	Mouse	
Sequence	<p>SQFNLSLEELGSNTGIQVFNQIIKSRPHENVVVSPHGIASI LGMLQLGA DGKTKKQLSTVMRYNVNGVGKVLKKINKAIVSKKNKDITV ANAVFLRNG FKMEVPFAVRNKDVFQCEVQNVNFQDPASASESINFWV KNETRGMDNLL SPNLIDGALTRLVLVNAVYFKGLWKSRLFQPESTKKRTFVA GDGKSYQVPM LAQLSVFRSGSTRTPNGLWYNFIELPYHGESISMLIALPTES STPLSAIL PHITTKTIDSWMNTMVPKRMQLVLPKFTAVAQTDLKEPLK ALGITEMFEP SKANFTKITRSESLHVSHILQKAKIEVSEDGKASAATTAILI ARSSPPW FVDRPFLFSIRHNPTGAILFLGQVKNKPLEHHHHHHH</p>	
Predicted molecular weight	43 kDa including tags	
Amino acids	20 to 397	
Tags	His tag C-Terminus	
Additional sequence information	This product is the mature full length protein from aa 20 to 397. The signal peptide is not included (NP_033281).	

Specifications

Our [Abpromise guarantee](#) covers the use of **ab222348** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE
Form	Liquid
Additional notes	Previously labelled as SERPINE2.


Preparation and Storage

Stability and Storage	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. pH: 7.40 Constituents: 10% Glycerol (glycerin, glycerine), PBS
------------------------------	---

General Info

Function	Serine protease inhibitor with activity toward thrombin, trypsin, and urokinase. Promotes neurite extension by inhibiting thrombin. Binds heparin.
Sequence similarities	Belongs to the serpin family.
Cellular localization	Secreted > extracellular space.

Images



15% SDS-PAGE analysis of 3 µg ab222348.

MWt 40-57 kDa (SDS-PAGE under reducing conditions).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors