

Product datasheet

Recombinant Outer capsid glycoprotein VP7 (Tagged)
 ab238289

1 Image

Description

Product name	Recombinant Outer capsid glycoprotein VP7 (Tagged)
Purity	> 85 % SDS-PAGE.
Expression system	Escherichia coli
Accession	A8D8S8
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Sequence	<pre> QNYGVNLPITGSMDTAYANSTQSEPFLTSTLCLYYPVEASN EIADEWKD TLSQLFLTKGWPTGSVYLKEYADIAAFSVEPQLYCDYNLVL MKYDSTQEL DMSELADLILNEWLCNPMIDITLYYQQTDEANKWISMGSS CTVKVCPLNT QTLGIGCLITNPDTFETVATTEKLVITDVVDGVSHKLNVTTA TCTIRNCK KLGPKENVAVIQVGGANILDITADPTTTPQTERMMAIWKK WWQVVYPVV DYVNQIIQTMSKRSLNSSAFYYRV </pre>
Predicted molecular weight	36 kDa including tags
Amino acids	34 to 309
Additional sequence information	N-terminal 10xHis-tagged and C-terminal Myc-tagged. Rotavirus A (strain RVA/Cow/Canada/C486/1977/G6P6[1]) (RV-A).

Specifications

Our [Abpromise guarantee](#) covers the use of **ab238289** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications SDS-PAGE

Form Liquid


Preparation and Storage

Stability and Storage

Shipped at 4°C. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

Constituents: Tris buffer, 50% Glycerol (glycerin, glycerine)

Images



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel analysis of ab238289.

SDS-PAGE - Recombinant Outer capsid glycoprotein VP7 (Tagged) (ab238289)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors