

Product datasheet

Recombinant Penicillin-binding protein 1A (His tag)
 ab225616

1 Image

Description

Product name	Recombinant Penicillin-binding protein 1A (His tag)	
Purity	> 90 % SDS-PAGE.	
Expression system	Escherichia coli	
Accession	A516G4	
Protein length	Protein fragment	
Animal free	No	
Nature	Recombinant	
Sequence	VDRISGKLPTQLSYRDPRGSTVYNEFFINGTIPTTEYDDIHVE AQINKLTG KLASKFTSPFLVESRVFLRRDYSPGVELLDQQWLLPYSID EGGSLPPTEE KNNSNTRDKNKDKNKNKNKDKNPSQDKPNNNNNDNNSN NNNNNNDNNNNT KPPENDSNQNHEDNKNKQ	
Predicted molecular weight	35 kDa including tags	
Amino acids	663 to 830	
Tags	His tag N-Terminus	
Additional sequence information	Clostridium botulinum (strain Hall / ATCC 3502 / NCTC 13319 / Type A) protein with 6xHis-SUMO tag at the N-terminus.	

Specifications

Our [Abpromise guarantee](#) covers the use of **ab225616** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE
Form	Liquid

Preparation and Storage

Stability and Storage Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.


Constituents: 50% Glycerol (glycerin, glycerine), Tris buffer

General Info

Cellular localization

Cell inner membrane

Images



ab225616 analyzed by SDS-PAGE (Tris-Glycine gel with 5% enrichment gel and 15% separation gel).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors