

Product datasheet

Recombinant Pig IL-8 protein ab90394

Description

Product name	Recombinant Pig IL-8 protein
Purity	> 95 % SDS-PAGE. Purified by Ion-exchange chromatography
Expression system	Yeast
Protein length	Full length protein
Animal free	No
Nature	Recombinant
Species	Pig
Sequence	ARVSAELRCQ CINTHSTPFH PKFIKELRVI ESGPHCENSE IMKLVNGKE VCLDPKEKWV QKVVQIFLKR TEKQQQQQ

Specifications

Our [Abpromise guarantee](#) covers the use of **ab90394** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE
Form	Lyophilised
Additional notes	Reconstituted protein is stable for at least 3 months when stored in working aliquots with a carrier protein at -20°C. Avoid repeated freeze/thaw cycles.

Preparation and Storage

Stability and Storage	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C. Constituents: 2.68% PBS, 10% Trehalose
Reconstitution	Reconstitute with sterile phosphate-buffered saline containing at least 0.1% carrier protein.

General Info

Function	IL-8 is a chemotactic factor that attracts neutrophils, basophils, and T-cells, but not monocytes. It is also involved in neutrophil activation. It is released from several cell types in response to an
-----------------	---

inflammatory stimulus. IL-8(6-77) has a 5-10-fold higher activity on neutrophil activation, IL-8(5-77) has increased activity on neutrophil activation and IL-8(7-77) has a higher affinity to receptors CXCR1 and CXCR2 as compared to IL-8(1-77), respectively.

Sequence similarities

Belongs to the intercrine alpha (chemokine CxC) family.

Post-translational modifications

Several N-terminal processed forms are produced by proteolytic cleavage after secretion from at least peripheral blood monocytes, leukocytes and endothelial cells. In general, IL-8(1-77) is referred to as interleukin-8. IL-8(6-77) is the most prominent form.

Cellular localization

Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors