

Product datasheet

Recombinant Protein ORF3 (Tagged) ab238356

1 Image

Description

Product name	Recombinant Protein ORF3 (Tagged)	
Purity	> 90 % SDS-PAGE.	
Expression system	Escherichia coli	
Accession	O90299	
Protein length	Full length protein	
Animal free	No	
Nature	Recombinant	
Species	Hepatitis E virus	
Sequence	MGSRPWALGLFCCCSSCFCLCCSRHRPVSRLAAVV GGAAAVPAVVSGVTG LILSPSQSPIFIQPTSPRMSPLRPGLDLVFANPSDHSA PLGATRPSAPP LPHVVDLPQLGPRR	
Predicted molecular weight	28 kDa including tags	
Amino acids	1 to 114	
Tags	His tag N-Terminus	
Additional sequence information	Genotype 1 (isolate Human/India/Hyderabad) (HEV-1). N-terminal 6xHis-SUMO-tag.	

Specifications

Our [Abpromise guarantee](#) covers the use of **ab238356** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	SDS-PAGE
Form	Liquid

Preparation and Storage


Stability and Storage	Shipped at 4°C. Store at -20°C or -80°C. Avoid freeze / thaw cycle. Constituents: Tris buffer, 50% Glycerol
------------------------------	--

General Info

Relevance

May act as a viral regulatory protein involved in the modulation of mitogenic signaling pathways. May be involved in virion morphogenesis and viral pathogenesis. Expedites the processing and secretion of AMBP from the hepatocyte.

Images



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) analysis with 5% enrichment gel and 15% separation gel of ab238356.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors