

Product datasheet

RS Repeat Peptide peptide ab204878

Description

Product name	RS Repeat Peptide peptide
Animal free	No
Nature	Synthetic
Sequence	GRSRSRSRSRSRSRSR

Specifications

Our [Abpromise guarantee](#) covers the use of **ab204878** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Applications	HPLC Functional Studies
Form	Lyophilised
Additional notes	ab204878 (RS Repeat Peptide peptide) can be utilized as a substrate for the following active protein kinases:

- [ab133149](#) (Active human BubR1 protein fragment)
- [ab131665](#) (Active human EIF2AK1 full length protein)

Preparation and Storage

Stability and Storage	Shipped at 4°C. Store at -20°C. Avoid freeze / thaw cycle.
Reconstitution	Dilute peptide in 20mM Tris-HCl, pH 7.5 solution to a final concentration of 1 mg/ml. For optimal storage, aliquot diluted product into smaller quantities and store at recommended temperature.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors