


Product datasheet

(R,S)-Sulforaphane, antagonist of Ah receptor ab141969

1 References 1 Image

Overview


<b>Product name</b>	(R,S)-Sulforaphane, antagonist of Ah receptor
<b>Description</b>	Potent, non-competitive antagonist of Ah receptor
<b>Biological description</b>	Potent, non-competitive antagonist of the aryl hydrocarbon (Ah) receptor (IC <sub>50</sub> = 1 μM). Isothiocyanate. Induces epoxide hydrolase. Anticarcinogenic. Active <i>in vivo</i> and <i>in vitro</i> .
<b>Purity</b>	> 98%
<b>General notes</b>	This compound is a neat liquid at room temperature, clear to slightly yellow in color. It is not dissolved in solvent. In small quantities the material will coat the sides of the ampule it is shipped in and may appear to be invisible, or it may be trapped in the tip of the ampule. Wash the upper and lower parts of the ampule with solvent to ensure all of the material is obtained.
<b>CAS Number</b>	4478-93-7
<b>Chemical structure</b>	

Properties

<b>Chemical name</b>	1-Isothiocyanato-4-methylsulfanylbutane
<b>Molecular weight</b>	177.29
<b>Molecular formula</b>	C <sub>6</sub> H <sub>11</sub> NOS <sub>2</sub>
<b>PubChem identifier</b>	5350
<b>Storage instructions</b>	Store at -20°C. Store under desiccating conditions. The product can be stored for up to 12 months.
<b>Solubility overview</b>	Soluble in DMSO to 100 mM
<b>Handling</b>	Wherever possible, you should prepare and use solutions on the same day. However, if you need to make up stock solutions in advance, we recommend that you store the solution as aliquots in tightly sealed vials at -20°C. Generally, these will be useable for up to one month. Before use, and prior to opening the vial we recommend that you allow your product to equilibrate to room temperature for at least 1 hour.  Need more advice on solubility, usage and handling? Please visit our <a href="#">frequently asked questions (FAQ) page</a> for more details.

<b>SMILES</b>	CS(=O)CCCCN=C=S
<b>Source</b>	Synthetic

## Images



2D chemical structure image of ab141969, (R,S)-Sulforaphane, antagonist of Ah receptor

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES, NOT FOR USE IN HUMANS"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors
- Abcam biochemicals are novel compounds and we have not tested their biological activity in house. Please use the literature to identify how to use these products effectively. If you require further assistance please contact the scientific support team