



Product datasheet

Anti-S tag antibody (Agarose) ab19369

Overview

Product name	Anti-S tag antibody (Agarose)
Description	Goat polyclonal to S tag (Agarose)
Host species	Goat
Conjugation	Agarose
Tested applications	Suitable for: IP
Species reactivity	Reacts with: Species independent
Immunogen	Synthetic peptide: KETAAAKFERQHMS (S tag) conjugated to KLH.
	 Run BLAST with  Run BLAST with
General notes	Affinity purified antibodies were coupled to agarose beads using a cyanogen bromide method. 100 µg antibody/200 µl gel. Volume: 400 µl of 50% slurry containing 200 µl gel.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C.
Storage buffer	pH: 7.20 Preservative: 0.1% Sodium azide Constituents: 0.0268% PBS, 0.58% Sodium chloride
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab19369** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IP		Use at an assay dependent dilution.

Target

Relevance S tag is an epitope tag composed of a 15 residue peptide, KETAAAKFERQHMS, derived from the pancreatic ribonuclease A.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors