

Product datasheet

Salmon Sperm ssDNA ab229278

Overview

Product name	Salmon Sperm ssDNA
General notes	<p>This product is a convenient, ready-to-use solution of salmon sperm DNA especially prepared for use in the preparation of pre-hybridization and hybridization solutions and as a DNA carrier in yeast transformation protocols and other related methods. This solution contains sheared single-stranded DNA molecules that can be used to block the non-specific attachment of probe DNA to the surface of a membrane (Southern) or to increase yeast transformation efficiency.</p> <p>Background: Prepared by a modification of the method of Emanuel and Chaikoff, JBC 203, 164 (1953) from salmon testes DNA by mechanical shearing and heat denaturation to an average size of 100 to 200 base pairs. To reverse any renaturation occurring during storage, this material should be briefly boiled and rapidly chilled prior to use. >75% native nucleic acid.</p> <p>Purity/Specificity: A260/A280: 1.85 % Hyperchromicity: 30.5.</p>
Tested applications	Suitable for: ChIP, In situ hybridization, Blocking

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
Storage buffer	Constituents: 0.16% Tris HCl, 0.03% EDTA
Relevance	The cell death process during apoptosis is essential for normal tissue homeostasis and cellular differentiation in multi cellular organisms. Single stranded DNA (ss DNA) is produced by nuclear fragmentation during the early stages of the apoptotic process.
Cellular localization	Nuclear

Applications

Our [Abpromise guarantee](#) covers the use of **ab229278** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ChIP		Use at an assay dependent concentration.

Application	Abreviews	Notes
In situ hybridization		Use at an assay dependent concentration.
Blocking		Use at an assay dependent concentration.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors