

Product datasheet

Anti-SAMD9L antibody ab200359

2 Images

Overview

Product name	Anti-SAMD9L antibody
Description	Rabbit polyclonal to SAMD9L
Host species	Rabbit
Tested applications	Suitable for: WB, IHC-P
Species reactivity	Reacts with: Human
Immunogen	Synthetic peptide corresponding to Human SAMD9L. (Near the N terminal; NP_689916). Database link: Q8IVG5
Positive control	Human colon cancer tissue lysate; Human liver cancer tissue.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
Storage buffer	pH: 7.3 Preservative: 0.05% Sodium azide Constituents: 50% Glycerol, 49% PBS
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab200359** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/200 - 1/1000. Predicted molecular weight: 185 kDa.
IHC-P		1/25 - 1/100. Perform heat mediated antigen retrieval before commencing with IHC staining protocol.

Target


Relevance

The function of this protein has not been specifically defined.

Cellular localization

early endosome


Images



Anti-SAMD9L antibody (ab200359) at 1/300 dilution + Human colon cancer tissue lysate at 40 µg

Predicted band size: 185 kDa

6% SDS-PAGE.



Immunohistochemical analysis of paraffin-embedded Human liver cancer tissue labeling SAMD9L with ab200359 at 1/35 dilution.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors