

Product datasheet

Anti-Smad2 + Smad3 antibody [EPR19557] - ChIP Grade ab207447

Recombinant RabMAb

[6 References](#) [10 Images](#)

Overview

Product name	Anti-Smad2 + Smad3 antibody [EPR19557] - ChIP Grade
Description	Rabbit monoclonal [EPR19557] to Smad2 + Smad3 - ChIP Grade
Host species	Rabbit
Tested applications	Suitable for: ICC/IF, IP, Flow Cyt, ChIP, WB
Species reactivity	Reacts with: Human
Immunogen	Recombinant fragment within Human Smad2 + Smad3 aa 1-300. The exact sequence is proprietary. Database link: Q15796
Positive control	WB: Recombinant protein fragment human Smad2; recombinant protein fragment human Smad3; HeLa, K562, A549 and Jurkat whole cell lysates; human fetal liver lysate. ICC/IF: HeLa cells. Flow Cyt: HeLa cells. IP: HeLa whole cell lysate. ChIP: Chromatin was prepared from HaCaT cells treated with 7ng/ml TGF- β for 1h.

General notes

This product is a recombinant monoclonal antibody, which offers several advantages including:

- High batch-to-batch consistency and reproducibility
- Improved sensitivity and specificity
- Long-term security of supply
- Animal-free production

For more information [see here](#).

Our RabMAb[®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to [RabMAb[®] patents](#).

Reproducibility is key to advancing scientific discovery and accelerating scientists' next breakthrough.

Abcam is leading the way with our range of recombinant antibodies, knockout-validated antibodies and knockout cell lines, all of which support improved reproducibility.

We are also planning to innovate the way in which we present recommended applications and species on our product datasheets, so that only applications & species that have been tested in our own labs, our suppliers or by selected trusted collaborators are covered by our Abpromise[™]

guarantee.

In preparation for this, we have started to update the applications & species that this product is Abpromise guaranteed for.

We are also updating the applications & species that this product has been “predicted to work with,” however this information is not covered by our Abpromise guarantee.

Applications & species from publications and Abreviews that have not been tested in our own labs or in those of our suppliers are not covered by the Abpromise guarantee.

Please check that this product meets your needs before purchasing. If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, as well as customer reviews and Q&As.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
Storage buffer	pH: 7.2 Preservative: 0.01% Sodium azide Constituents: PBS, 40% Glycerol (glycerin, glycerine), 0.05% BSA
Purity	Protein A purified
Clonality	Monoclonal
Clone number	EPR19557
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab207447** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ICC/IF		1/500.
IP		1/40.
Flow Cyt		1/500.
ChIP		Use 2 µg for 25 µg of chromatin.
WB		1/2000. Detects a band of approximately 58-62 kDa (predicted molecular weight: 52 kDa).

Target


Relevance SMAD is a family of proteins similar to the gene products of the Drosophila gene 'mothers against

decapentaplegic' (Mad) and the *C. elegans* gene Sma. SMAD proteins are signal transducers and transcriptional modulators that mediate multiple signaling pathways. They mediate the signal of the transforming growth factor (TGF)-beta, and thus regulate multiple cellular processes, such as cell proliferation, apoptosis, and differentiation.

Cellular localization

Cytoplasm. Nucleus. Note: Cytoplasmic in the absence of ligand. Migrates to the nucleus when complexed with SMAD4.

Images



All lanes : Anti-Smad2 + Smad3 antibody [EPR19557] - ChIP Grade (ab207447) at 1/50000 dilution

Lane 1 : Recombinant protein fragment human Smad2

Lane 2 : Recombinant protein fragment human Smad3

Lysates/proteins at 0.01 µg per lane.

Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) (ab97051) at 1/100000 dilution

Predicted band size: 52 kDa


Observed band size: 38,60 kDa

[why is the actual band size different from the predicted?](#)

Blocking/Dilution buffer: 5% NFDm/TBST.

Exposure time: Lane 1: 2 seconds; Lane 2: 15 seconds.

Human Smad2 fragment recombinant protein contains aa2-270 with His-Tag[®]. Human Smad3 fragment recombinant protein contains aa2-227 with His-Tag[®] and GST-tag.



Western blot - Anti-Smad2 + Smad3 antibody [EPR19557] - ChIP Grade (ab207447)

All lanes : Anti-Smad2 + Smad3 antibody [EPR19557] - ChIP Grade (ab207447) at 1/10000 dilution

Lane 1 : HeLa (Human T cell leukemia cell line from peripheral blood) whole cell lysate

Lane 2 : K562 (Human chronic myelogenous leukemia cell line from bone marrow) whole cell lysate

Lane 3 : A549 (Human lung carcinoma cell line) whole cell lysate

Lane 4 : Human fetal liver lysate

Lysates/proteins at 10 µg per lane.

Secondary


All lanes : Goat Anti-Rabbit IgG H&L (HRP) (ab97051) at 1/100000 dilution

Predicted band size: 52 kDa

Observed band size: 58-62 kDa [why is the actual band size different from the predicted?](#)

Exposure time: 3 minutes

Blocking/Dilution buffer: 5% NFDm/TBST.



Western blot - Anti-Smad2 + Smad3 antibody [EPR19557] - ChIP Grade (ab207447)

Anti-Smad2 + Smad3 antibody [EPR19557] - ChIP Grade (ab207447) at 1/2000 dilution + Jurkat (Human T cell leukemia cell line from peripheral blood) whole cell lysate at 10 µg

Secondary


Goat Anti-Rabbit IgG H&L (HRP) (ab97051) at 1/100000 dilution

Predicted band size: 52 kDa

Observed band size: 58-62 kDa [why is the actual band size different from the predicted?](#)

Exposure time: 15 seconds

Blocking/Dilution buffer: 5% NFDm/TBST.



Western blot - Anti-Smad2 + Smad3 antibody [EPR19557] - ChIP Grade (ab207447)

Anti-Smad2 + Smad3 antibody [EPR19557] - ChIP Grade (ab207447) at 1/10000 dilution + HeLa (Human epithelial cell line from cervix adenocarcinoma) whole cell lysate at 10 μ g

Secondary

Goat Anti-Rabbit IgG H&L (HRP) (ab97051) at 1/100000 dilution


Predicted band size: 52 kDa

Observed band size: 52,58-62 kDa [why is the actual band size different from the predicted?](#)

Exposure time: 30 seconds

Blocking/Dilution buffer: 5% NFDM/TBST.

Smad2 and Smad3 can be resolved using a lower percentage gel.



Immunocytochemistry/ Immunofluorescence - Anti-Smad2 + Smad3 antibody [EPR19557] - ChIP Grade (ab207447)


Immunofluorescent analysis of 4% paraformaldehyde-fixed, 0.1% Triton X-100 permeabilized HeLa (Human epithelial cell line from cervix adenocarcinoma) cells labeling Smad2 + Smad3 with ab207447 at 1/500 dilution, followed by Goat Anti-Rabbit IgG (Alexa Fluor[®] 488) (ab150077) secondary antibody at 1/1000 dilution (green). Confocal image showing signal translocation from cytoplasm to nucleus after TGF-beta (10ng/ml, 1h) treatment in HeLa cells. PMID: 9006934. The nuclear counter stain is DAPI (blue).

Tubulin is detected with Anti-alpha Tubulin mouse MAb (ab7291) at 1/1000 dilution, followed by Goat Anti-Mouse IgG H&L (Alexa Fluor[®] 594) (ab150120) secondary antibody at 1/1000 dilution (red).

The negative controls are as follows:


-ve control 1: ab207447 at 1/500 dilution, followed by Goat Anti-Mouse IgG H&L (Alexa Fluor[®] 594) (ab150120) secondary at 1/1000 dilution.

-ve control 2: Anti-alpha Tubulin mouse MAb (ab7291) at 1/1000 dilution, followed by Goat Anti-Rabbit IgG H&L (Alexa Fluor[®] 488) (ab150077) secondary at 1/1000 dilution.



Flow Cytometry - Anti-Smad2 + Smad3 antibody
[EPR19557] - ChIP Grade (ab207447)

Flow cytometric analysis of 4% paraformaldehyde-fixed HeLa (Human epithelial cell line from cervix adenocarcinoma) cells labeling Smad2 and Smad3 with ab207447 at 1/500 dilution (red) compared with a rabbit monoclonal IgG isotype control (ab172730; black) and an unlabelled control (cells without incubation with primary antibody and secondary antibody; blue). Goat anti rabbit IgG (Alexa Fluor® 488) at 1/500 dilution was used as the secondary antibody.



Immunoprecipitation - Anti-Smad2 + Smad3
antibody [EPR19557] - ChIP Grade (ab207447)

Smad2 + Smad3 were immunoprecipitated from 1mg of HeLa (Human epithelial cell line from cervix adenocarcinoma) whole cell lysate with ab207447 at 1/40 dilution. Western blot was performed from the immunoprecipitate using ab207447 at 1/1000 dilution. VeriBlot for IP Detection Reagent (HRP) (ab131366), was used for detection at 1/10000 dilution.


Lane 1: HeLa whole cell lysate 10µg (Input).

Lane 2: ab207447 IP in HeLa whole cell lysate.

Lane 3: Rabbit monoclonal IgG (ab172730) instead of ab207447 in HeLa whole cell lysate.

Blocking and dilution buffer and concentration: 5% NFDM/TBST.

Exposure time: 5 seconds.



Immunoprecipitation - Anti-Smad2 + Smad3 antibody [EPR19557] - ChIP Grade (ab207447)

ab207447 Immunoprecipitating Smad2 + Smad3 in human HeLa whole cell lysate . 2µg of capture antibody in 0.35mg lysate was used. 10µg of cell lysate was incubated with primary antibody at a dilution of 1/1000 and VeriBlot for IP Detection Reagent (HRP) [ab131366](#) was used for detection at a dilution of 1/1000.

Lane 1: HeLa (human cervix adenocarcinoma) whole cell lysate


Lane 2: HeLa (human cervix adenocarcinoma) whole cell lysate

Lane 3: Rabbit monoclonal IgG ([ab172730](#)) instead of ab207447 in HeLa (human cervix adenocarcinoma) whole cell lysate

Blocking buffer: 5% NFDm/TBST

Diluting buffer: 5% NFDm/TBST

Smad2 and Smad3 can be resolved using a lower percentage gel.



Every new batch of this antibody is tested at Abcam in ChIP.

ChIP - Anti-Smad2 + Smad3 antibody [EPR19557] - ChIP Grade (ab207447)

Chromatin was prepared from HaCaT (Human keratinocyte cell line) cells treated with 7ng/ml TGF-β for 1h and non-treated according to the Abcam X-ChIP protocol. Cells were fixed with formaldehyde for 10 minutes. The ChIP was performed with 25µg of chromatin, 2µg of ab207447 (blue), and 20µl of Anti rabbit IgG sepharose beads. 2µg of rabbit normal IgG was added to the beads control (yellow). The immunoprecipitated DNA was quantified by real time PCR (Sybr green approach).

The ChIP condition is designed against Smad2 refer to PMID: 18955504.

Why choose a recombinant antibody?



Research with confidence
Consistent and reproducible results



Long-term and scalable supply
Recombinant technology



Success from the first experiment
Confirmed specificity



Ethical standards compliant
Animal-free production

Anti-Smad2 + Smad3 antibody [EPR19557] - ChIP Grade (ab207447)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors