




Product datasheet

Anti-SNRPB2 antibody ab234872

1 Image

Overview

<b>Product name</b>	Anti-SNRPB2 antibody
<b>Description</b>	Rabbit polyclonal to SNRPB2
<b>Host species</b>	Rabbit
<b>Tested applications</b>	<b>Suitable for:</b> IHC-P
<b>Species reactivity</b>	<b>Reacts with:</b> Human <b>Predicted to work with:</b> Mouse 
<b>Immunogen</b>	Recombinant fragment corresponding to Human SNRPB2 aa 1-150. Sequence: MDIRPNHTYINMNDKIKKEELKRSLYALFSQFGHVVDI VALKTMKMRG QAFVIFKELGSSTNALRQLQGFPFYGKPMRIQYAKTDS DISKMRGTFAD KEKKKEKKKAKTVEQTATTTNKKPGQGTPNSANTQGN STPNPQVPDYPPN  Database link: <a href="#">SNRPB2</a>   <a href="#">Run BLAST with</a>  <a href="#">Run BLAST with</a>
<b>Positive control</b>	IHC-P: Human cervical cancer tissue.

Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
<b>Storage buffer</b>	pH: 7.30 Preservative: 0.02% Sodium azide Constituents: PBS, 50% Glycerol
<b>Purity</b>	Immunogen affinity purified
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG

## Applications

Our [Abpromise guarantee](#) covers the use of **ab234872** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IHC-P		1/20 - 1/200.

## Target

<b>Function</b>	Involved in pre-mRNA splicing. This protein is associated with snRNP U2. It binds stem loop IV of U2 snRNA only in presence of the U2A' protein.
<b>Sequence similarities</b>	Belongs to the RRM U1 A/B" family. Contains 2 RRM (RNA recognition motif) domains.
<b>Cellular localization</b>	Nucleus.

## Images



Paraffin-embedded human cervical cancer tissue stained for SNRNP2 using ab234872 at 1/100 dilution in immunohistochemical analysis.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-SNRNP2 antibody (ab234872)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors