




Product datasheet

Anti-SNX7 antibody ab220662

1 Image

Overview

Product name	Anti-SNX7 antibody
Description	Rabbit polyclonal to SNX7
Host species	Rabbit
Tested applications	Suitable for: ICC/IF
Species reactivity	Reacts with: Human Predicted to work with: Mouse, Rat, Cynomolgus monkey 
Immunogen	Recombinant fragment corresponding to Human SNX7 aa 227-301. Sequence: IHILWSASEEDLVDTLKDVASCIDRCCKATEKRMSGLS EALLPVVHEYVL YSEMLMGVMKRRDQIQAEELDSKVEV Database link: Q9UNH6  Run BLAST with  Run BLAST with
Positive control	A431 cells.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
Storage buffer	pH: 7.2 Preservative: 0.02% Sodium azide Constituents: 40% Glycerol, 59% PBS
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab220662** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ICC/IF		Use a concentration of 1 - 4 µg/ml. Use PFA/Triton X-100 Fixation/Permeabilization.

Target

Function	May be involved in several stages of intracellular trafficking.
Sequence similarities	Belongs to the sorting nexin family. Contains 1 PX (phox homology) domain.
Cellular localization	Cytoplasmic vesicle membrane.

Images



Immunofluorescent analysis of A431 cells (PFA-fixed/Triton X-100 permeabilized) labeling SNX7 with ab220662 at 4µg/ml (green).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors