

Product datasheet

Streptavidin-HRP Solution ab210901

1 References

Overview

Product name Streptavidin-HRP Solution

Product overview Streptavidin-HRP Solution is recommended for use with the Matched Antibody Pair kits for detecting biotinylated molecules during ELISA development.

See our [website](#) for a protocol.

We recommend an initial 1:500 - 1:5,000 dilution range of the Streptavidin-HRP Solution. Optimal dilution should however be determined by the end-user.

The Streptavidin-HRP Solution is also available as part of our Matched Antibody Pair Kit Accessory Pack ([ab210905](#)).

Concentration: 50 µg/mL

Properties

Storage instructions Store at -20°C. Please refer to protocols.

Components	110 µl
Streptavidin HRP Solution (ab210901) 50 µg/mL	1 x 110µl

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors