

## Product datasheet

# TEO-Tricine Precast Gels - RunBlue™ (4-20%, 2-well, 10x10cm) ab270461

3 Images

### Overview

<b>Product name</b>	TEO-Tricine Precast Gels - RunBlue™ (4-20%, 2-well, 10x10cm)
<b>Tested applications</b>	<b>Suitable for:</b> SDS-PAGE
<b>General notes</b>	<p>TEO-Tricine Precast Gels - RunBlue™ (4-20%, 2-well, 10x10cm) offers superior rigidity and stability over traditional polyacrylamide gels and are compatible with the most common tanks and systems. For convenience the precast gels come with the comb already removed as the cassette locks the fingers in place, for additional convenience there is no tape to be removed and cassettes can be opened by hand.</p> <p>10 gel pack.</p> <p>This product is manufactured by Expedeon, an Abcam company. It was previously called RunBlue™ TEO-Tricine SDS Gel 4-20% 2 well - 10x10. Expedeon product code NXG42002 is the same as the 1 pack size of this product.</p>


### Properties

<b>Form</b>	Solid
<b>Storage instructions</b>	Shipped at Room Temperature. Store at Room Temperature. Store at +4°C. Do Not Freeze.
<b>Storage buffer</b>	<p>pH: 7.8</p> <p>Constituents: 0.2% Sodium lauryl sulfate, 5% Triethanolamine, 20% Polyacrylamide, 0.02% Acrylamide, 0.02% Bis-acrylamide</p>


### Applications

**The Abpromise guarantee** Our [Abpromise guarantee](#) covers the use of ab270461 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
SDS-PAGE		Use at an assay dependent concentration.




SDS-PAGE - TEO-Tricine Precast Gels - RunBlue™  
(4-20%, 2-well, 10x10cm) (ab270461)



SDS-PAGE - TEO-Tricine Precast Gels - RunBlue™  
(4-20%, 2-well, 10x10cm) (ab270461)

## Gel benefits - Click to expand



SDS-PAGE - TEO-Tricine Precast Gels - RunBlue™  
(4-20%, 2-well, 10x10cm) (ab270461)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors