

## Product datasheet

# Anti-TGF beta Receptor II antibody [EPR23148-21] ab269279

Recombinant RabMAb

★★★★★ 1 Abreviews 1 References 4 Images

### Overview

<b>Product name</b>	Anti-TGF beta Receptor II antibody [EPR23148-21]
<b>Description</b>	Rabbit monoclonal [EPR23148-21] to TGF beta Receptor II
<b>Host species</b>	Rabbit
<b>Tested applications</b>	<b>Suitable for:</b> WB, IP <b>Unsuitable for:</b> Flow Cyt, ICC/IF or IHC-P
<b>Species reactivity</b>	<b>Reacts with:</b> Mouse
<b>Immunogen</b>	Recombinant fragment. This information is proprietary to Abcam and/or its suppliers.
<b>Positive control</b>	WB: Mouse lung, Mouse spleen, Hepa1-6, bEnd.3, EL4, IL-2, Mouse thymus and Mouse lung, Mouse lung treated with PNGase F tissue lysates. IP: Mouse lung tissue lysate.
<b>General notes</b>	This product is a recombinant monoclonal antibody, which offers several advantages including: <ul style="list-style-type: none"> <li>- High batch-to-batch consistency and reproducibility</li> <li>- Improved sensitivity and specificity</li> <li>- Long-term security of supply</li> <li>- Animal-free production</li> </ul> For more information <a href="#">see here</a> . Our RabMAb <sup>®</sup> technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to <a href="#">RabMAb<sup>®</sup> patents</a> .

### Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
<b>Storage buffer</b>	pH: 7.2 Preservative: 0.01% Sodium azide Constituents: PBS, 40% Glycerol (glycerin, glycerine), 0.05% BSA
<b>Purity</b>	Protein A purified
<b>Clonality</b>	Monoclonal

Clone number                      EPR23148-21

Isotype                                IgG

## Applications

**The Abpromise guarantee**            Our [Abpromise guarantee](#) covers the use of ab269279 in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB	★★★★★ (1)	1/1000. Predicted molecular weight: 65 kDa.
IP		1/30.

**Application notes**                      Is unsuitable for Flow Cyt, ICC/IF or IHC-P.

## Target

### Function

Transmembrane serine/threonine kinase forming with the TGF-beta type I serine/threonine kinase receptor, TGFBR1, the non-promiscuous receptor for the TGF-beta cytokines TGFB1, TGFB2 and TGFB3. Transduces the TGFB1, TGFB2 and TGFB3 signal from the cell surface to the cytoplasm and is thus regulating a plethora of physiological and pathological processes including cell cycle arrest in epithelial and hematopoietic cells, control of mesenchymal cell proliferation and differentiation, wound healing, extracellular matrix production, immunosuppression and carcinogenesis. The formation of the receptor complex composed of 2 TGFBR1 and 2 TGFBR2 molecules symmetrically bound to the cytokine dimer results in the phosphorylation and the activation of TGFBR1 by the constitutively active TGFBR2. Activated TGFBR1 phosphorylates SMAD2 which dissociates from the receptor and interacts with SMAD4. The SMAD2-SMAD4 complex is subsequently translocated to the nucleus where it modulates the transcription of the TGF-beta-regulated genes. This constitutes the canonical SMAD-dependent TGF-beta signaling cascade. Also involved in non-canonical, SMAD-independent TGF-beta signaling pathways.

### Involvement in disease

Defects in TGFBR2 are the cause of hereditary non-polyposis colorectal cancer type 6 (HNPCC6) [MIM:614331]. Mutations in more than one gene locus can be involved alone or in combination in the production of the HNPCC phenotype (also called Lynch syndrome). Most families with clinically recognized HNPCC have mutations in either MLH1 or MSH2 genes. HNPCC is an autosomal, dominantly inherited disease associated with marked increase in cancer susceptibility. It is characterized by a familial predisposition to early onset colorectal carcinoma (CRC) and extra-colonic cancers of the gastrointestinal, urological and female reproductive tracts. HNPCC is reported to be the most common form of inherited colorectal cancer in the Western world, and accounts for 15% of all colon cancers. Cancers in HNPCC originate within benign neoplastic polyps termed adenomas. Clinically, HNPCC is often divided into two subgroups. Type I: hereditary predisposition to colorectal cancer, a young age of onset, and carcinoma observed in the proximal colon. Type II: patients have an increased risk for cancers in certain tissues such as the uterus, ovary, breast, stomach, small intestine, skin, and larynx in addition to the colon. Diagnosis of classical HNPCC is based on the Amsterdam criteria: 3 or more relatives affected by colorectal cancer, one a first degree relative of the other two; 2 or more generation affected; 1 or more colorectal cancers presenting before 50 years of age; exclusion of hereditary polyposis syndromes. The term "suspected HNPCC" or "incomplete HNPCC" can be used to describe families who do not or only partially fulfill the Amsterdam

criteria, but in whom a genetic basis for colon cancer is strongly suspected. HNPCC6 is a type of colorectal cancer complying with the clinical criteria of HNPCC, except that the onset of cancer was beyond 50 years of age in all cases.

Defects in TGFBR2 are a cause of esophageal cancer (ESCR) [MIM:133239].

Defects in TGFBR2 are the cause of Loeys-Dietz syndrome type 1B (LDS1B) [MIM:610168]. LDS1 is an aortic aneurysm syndrome with widespread systemic involvement. The disorder is characterized by arterial tortuosity and aneurysms, craniosynostosis, hypertelorism, and bifid uvula or cleft palate. Other findings include exotropia, micrognathia and retrognathia, structural brain abnormalities, intellectual deficit, congenital heart disease, translucent skin, joint hyperlaxity and aneurysm with dissection throughout the arterial tree.

Defects in TGFBR2 are the cause of Loeys-Dietz syndrome type 2B (LDS2B) [MIM:610380]. An aortic aneurysm syndrome with widespread systemic involvement. Physical findings include prominent joint laxity, easy bruising, wide and atrophic scars, velvety and translucent skin with easily visible veins, spontaneous rupture of the spleen or bowel, diffuse arterial aneurysms and dissections, and catastrophic complications of pregnancy, including rupture of the gravid uterus and the arteries, either during pregnancy or in the immediate postpartum period. LDS2 is characterized by the absence of craniofacial abnormalities with the exception of bifid uvula that can be present in some patients. Note=TGFBR2 mutations Cys-460 and His-460 have been reported to be associated with thoracic aortic aneurysms and dissection (TAAD). This phenotype, also known as thoracic aortic aneurysms type 3 (AAT3), is distinguished from LDS2B by having aneurysms restricted to thoracic aorta. As individuals carrying these mutations also exhibit descending aortic disease and aneurysms of other arteries (PubMed:16027248), they have been considered as LDS2B by the OMIM resource.

### Sequence similarities

Belongs to the protein kinase superfamily. TKL Ser/Thr protein kinase family. TGFBR2 receptor subfamily.

Contains 1 protein kinase domain.


### Post-translational modifications

Phosphorylated on a Ser/Thr residue in the cytoplasmic domain.

### Cellular localization

Cell membrane.

## Images



Western blot - Anti-TGF beta Receptor II antibody  
[EPR23148-21] (ab269279)

**All lanes** : Anti-TGF beta Receptor II antibody [EPR23148-21]  
(ab269279) at 1/1000 dilution

**Lane 1** : Mouse lung tissue lysate

**Lane 2** : Mouse lung tissue lysate treated with PNGase F

Lysates/proteins at 20 µg per lane.

### Secondary

**All lanes** : Goat Anti-Rabbit IgG, (H+L), Peroxidase conjugated  
(ab97051) at 1/100000 dilution

**Predicted band size:** 65 kDa


**Observed band size:** 80/70/67 kDa

Blocking and diluting buffer and concentration: 5% NFDm/TBST

The expression profile observed is consistent with what has been described in the literature (PMID: 25323858).

Two bands of 80 and 70kDa are glycosylated forms. The band of 67kDa is de-glycosylated form.

Exposure time: 3 minutes



Immunoprecipitation - Anti-TGF beta Receptor II antibody [EPR23148-21] (ab269279)

TGF beta Receptor II was immunoprecipitated from 0.35 mg Mouse lung tissue lysate with ab269279 at 1/30 dilution (2ug in 0.35mg lysates). Western blot was performed on the immunoprecipitate using ab269279 at 1/500 dilution. VeriBlot for IP Detection Reagent (HRP)([ab131366](#)) was used at 1/1000 dilution.


Lane 1: Mouse lung tissue lysate 10 ug

Lane 2: ab269279 IP in Mouse lung tissue lysate

Lane 3: Rabbit monoclonal IgG ([ab172730](#)) instead of ab269279 in mouse lung tissue lysate

Blocking and dilution buffer and concentration: 5% NFDm/TBST.

Exposure time: 32 seconds



Western blot - Anti-TGF beta Receptor II antibody [EPR23148-21] (ab269279)

**All lanes :** Anti-TGF beta Receptor II antibody [EPR23148-21] (ab269279) at 1/1000 dilution

**Lane 1 :** Mouse lung tissue lysate

**Lane 2 :** Mouse spleen tissue lysate

**Lane 3 :** Hepa1-6 (mouse hepatoma epithelial cell) whole cell lysate

**Lane 4 :** bEnd.3 (mouse brain endothelioma ) whole cell lysate

**Lane 5 :** EL4.IL-2 (mouse lymphoma t lymphocyte) whole cell lysate

**Lane 6 :** Mouse thymus tissue lysate

Lysates/proteins at 40 µg per lane.

### Secondary

**All lanes :** Goat Anti-Rabbit IgG, (H+L), Peroxidase conjugated ([ab97051](#)) at 1/100000 dilution

**Predicted band size:** 65 kDa

**Observed band size:** 80/70 kDa

Blocking and diluting buffer and concentration: 5% NFD/MTBST

The expression profile observed is consistent with what has been described in the literature (PMID: 25323858).

Exposure time: Lanes 1-4: 70 seconds; Lanes 5-6: 3 minutes.

Why choose a recombinant antibody?

**Research with confidence**  
Consistent and reproducible results

**Long-term and scalable supply**  
Recombinant technology

**Success from the first experiment**  
Confirmed specificity

**Ethical standards compliant**  
Animal-free production

Anti-TGF beta Receptor II antibody [EPR23148-21]  
(ab269279)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors