

Product datasheet

Anti-TMIGD2/IGPR1 antibody ab121333

2 Images

Overview

<b>Product name</b>	Anti-TMIGD2/IGPR1 antibody
<b>Description</b>	Rabbit polyclonal to TMIGD2/IGPR1
<b>Host species</b>	Rabbit
<b>Tested applications</b>	<b>Suitable for:</b> IHC-P, WB
<b>Species reactivity</b>	<b>Reacts with:</b> Human
<b>Immunogen</b>	Recombinant fragment corresponding to Human TMIGD2/IGPR1 aa 44-146.
<b>Positive control</b>	Human duodenum tissue; TMIGD2/IGPR1 over-expression HEK293T lysate.
<b>General notes</b>	This product was previously labelled as TMIGD2

Properties

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
<b>Storage buffer</b>	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol
<b>Purity</b>	Immunogen affinity purified
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab121333** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IHC-P		1/50 - 1/200. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.


Application	Abreviews	Notes
WB		1/250 - 1/500.

## Target

**Sequence similarities** Contains 1 Ig-like (immunoglobulin-like) domain.

**Cellular localization** Membrane.

## Images




Western blot - Anti-TMIGD2/IGPR1 antibody (ab121333)

**All lanes** : Anti-TMIGD2/IGPR1 antibody (ab121333) at 1/250 dilution

**Lane 1** : Negative control (vector only transfected HEK293T lysate)

**Lane 2** : Over-expression Lysate (Co-expressed with a C-terminal myc-DDK tag (~3.1 kDa) in mammalian HEK293T cells)



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-TMIGD2/IGPR1 antibody (ab121333)

ab121333, at 1/50, staining TMIGD2/IGPR1 in paraffin embedded Human duodenum tissue by Immunohistochemistry.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors