

Product datasheet

Anti-TOX1 antibody [EP12540] - C-terminal ab181197

RabMAb

1 Image

Overview

Product name	Anti-TOX1 antibody [EP12540] - C-terminal
Description	Rabbit monoclonal [EP12540] to TOX1 - C-terminal
Host species	Rabbit
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Mouse, Rat, Human
Immunogen	Recombinant fragment within Human TOX1 aa 350 to the C-terminus. The exact sequence is proprietary. Database link: O94900
Positive control	Ramos, Jurkat and HuT-78 cell lysate.
General notes	Our RabMAb [®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb[®] patents .

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
Storage buffer	Preservative: 0.01% Sodium azide Constituents: 59% PBS, 40% Glycerol, 0.05% BSA
Purity	Tissue culture supernatant
Clonality	Monoclonal
Clone number	EP12540
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab181197** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/1000. Detects a band of approximately 58 kDa (predicted molecular weight: 58 kDa).

Target

Function	May play a role in regulating T-cell development.
Sequence similarities	Contains 1 HMG box DNA-binding domain.
Cellular localization	Nucleus.

Images



All lanes : Anti-TOX1 antibody [EP12540] - C-terminal (ab181197)
at 1/1000 dilution

Lane 1 : Ramos cell lysate

Lane 2 : Jurkat cell lysate

Lane 3 : Hut-78 cell lysate

Lysates/proteins at 10 µg per lane.

Secondary

All lanes : Goat Anti-Rabbit IgG, (H+L), Peroxidase conjugated at
1/1000 dilution

Predicted band size: 58 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors