

Product datasheet

# Anti-TRAF6BP/TAX1BP1 antibody [EPR13287(B)] - BSA and Azide free ab238935

Recombinant RabMAb

1 Image

Overview

<b>Product name</b>	Anti-TRAF6BP/TAX1BP1 antibody [EPR13287(B)] - BSA and Azide free
<b>Description</b>	Rabbit monoclonal [EPR13287(B)] to TRAF6BP/TAX1BP1 - BSA and Azide free
<b>Host species</b>	Rabbit
<b>Tested applications</b>	<b>Suitable for:</b> WB, ICC/IF <b>Unsuitable for:</b> Flow Cyt, IHC-P or IP
<b>Species reactivity</b>	<b>Reacts with:</b> Mouse, Rat, Human
<b>Immunogen</b>	Synthetic peptide within Human TRAF6BP/TAX1BP1 aa 750 to the C-terminus (C terminal) (Cysteine residue). The exact sequence is proprietary. Isoform 1 Database link: <a href="#">Q86VP1</a>
<b>Positive control</b>	ICC/IF: A549 cells.
<b>General notes</b>	ab238935 is the carrier-free version of <a href="#">ab176572</a> This format is designed for use in antibody labeling, including fluorochromes, metal isotopes, oligonucleotides, enzymes.

Our [carrier-free formats](#) are supplied in a buffer free of BSA, sodium azide and glycerol for higher conjugation efficiency.

Use our [conjugation kits](#) for antibody conjugates that are ready-to-use in as little as 20 minutes with <1 minute hands-on-time and 100% antibody recovery: available for fluorescent dyes, HRP, biotin and gold.

Ab238935 is compatible with the Maxpar® Antibody Labeling Kit from Fluidigm.

*Maxpar® is a trademark of Fluidigm Canada Inc.*

This product was previously labelled as TRAF6BP

This product is a recombinant monoclonal antibody, which offers several advantages including:

- High batch-to-batch consistency and reproducibility
- Improved sensitivity and specificity
- Long-term security of supply
- Animal-free production

For more information [see here](#).

Our RabMAb® technology is a patented hybridoma-based technology for making rabbit

monoclonal antibodies. For details on our patents, please refer to [RabMAb® patents](#).

## Properties

---

<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
<b>Storage buffer</b>	Constituent: PBS
<b>Carrier free</b>	Yes
<b>Purity</b>	Protein A purified
<b>Clonality</b>	Monoclonal
<b>Clone number</b>	EPR13287(B)
<b>Isotype</b>	IgG

## Applications

---

Our [Abpromise guarantee](#) covers the use of **ab238935** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		Use at an assay dependent concentration. Predicted molecular weight: 91 kDa.
ICC/IF		Use at an assay dependent concentration.

**Application notes**                      Is unsuitable for Flow Cyt, IHC-P or IP.

## Target

---


**Function**                                      Inhibits TNF-induced apoptosis by mediating the TNFAIP3 anti-apoptotic activity. Degraded by caspase-3-like family proteins upon TNF-induced apoptosis. May also play a role in the pro-inflammatory cytokine IL-1 signaling cascade.

**Tissue specificity**                      Expressed in all tissues tested.

---

## Images

---



Immunocytochemistry/ Immunofluorescence - Anti-TRAF6BP/TAX1BP1 antibody [EPR13287(B)] - BSA and Azide free (ab238935)

Immunofluorescence analysis of A549 (Human lung carcinoma cell line) cells labeling TRAF6BP/TAX1BP1 using [ab176572](#) at a 1/50 dilution (green). DAPI nuclear staining (blue).

This data was developed using the same antibody clone in a different buffer formulation containing PBS, BSA, glycerol, and sodium azide ([ab176572](#)).

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors