


Product datasheet

Anti-TRIOBP antibody - C-terminal ab151320

1 References 4 Images

Overview

Product name	Anti-TRIOBP antibody - C-terminal
Description	Rabbit polyclonal to TRIOBP - C-terminal
Host species	Rabbit
Tested applications	Suitable for: WB, IHC-P, ICC/IF
Species reactivity	Reacts with: Mouse, Human Predicted to work with: Cow 
Immunogen	Recombinant fragment, corresponding to a region within amino acids 1833 and 2061 of Human TRIOBP (Uniprot: Q9H2D6).
Positive control	A549, HeLa, HepG2, HCT116 and NIH 3T3 whole cell lysates; Human MDA MB 157 xenograft tissue; HeLa cells.
General notes	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
Storage buffer	pH: 7.00 Preservative: 0.01% Thimerosal (merthiolate) Constituents: 1.21% Tris, 0.75% Glycine, 20% Glycerol (glycerin, glycerine)
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgG

Applications

The Abpromise guarantee

Our [Abpromise guarantee](#) covers the use of ab151320 in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/500 - 1/3000. Predicted molecular weight: 261 kDa.
IHC-P		1/100 - 1/1000. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol. 10mM Citrate buffer
ICC/IF		1/100 - 1/1000.

Target


Relevance

TRIOBP is a protein that interacts with trio, which is involved with neural tissue development and controlling actin cytoskeleton organization, cell motility and cell growth. This trio binding protein also associates with F actin and stabilizes F actin structures.

Cellular localization

Nucleus. Cytoplasm; cytoskeleton. Note: Localized to F-actin in a periodic pattern.

Images



All lanes : Anti-TRIOBP antibody - C-terminal (ab151320) at 1/1000 dilution


Lane 1 : HeLa whole cell lysate

Lane 2 : HepG2 whole cell lysate

Lysates/proteins at 30 µg per lane.

Predicted band size: 261 kDa

10% SDS PAGE




Western blot - Anti-TRIOBP antibody - C-terminal (ab151320)

Anti-TRIOBP antibody - C-terminal (ab151320) at 1/1000 dilution + NIH 3T3 whole cell lysate at 30 μ g


Predicted band size: 261 kDa

10% SDS PAGE



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-TRIOBP antibody - C-terminal (ab151320)

Immunohistochemical analysis of paraffin-embedded Human MDA-MB-157 xenograft tissue labeling TRIOBP with ab151320 at 1/100 dilution.



Immunocytochemistry/ Immunofluorescence - Anti-TRIOBP antibody - C-terminal (ab151320)

Immunofluorescent analysis of paraformaldehyde-fixed HeLa cells labeling TRIOBP using ab151320 at 1/200 dilution. Right image shows cells co-stained with Hoechst 33342.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors