

Product datasheet

Tween 20 (Polyoxyethelenesorbitan Monolaurate) ab128987

1 References

Overview

Product name	Tween 20 (Polyoxyethelenesorbitan Monolaurate)
Tested applications	Suitable for: IHC-P
General notes	Tween 20 detergent (surfactant) ab128987 formulated to be added to buffers and reagents for immunohistochemistry to decrease background staining and enhance reagent spreading during IHC in automated and manual procedures. Procedure: Suggested to use at concentrations between 0.05% and 0.1% in buffers and reagents. IHC Protocol Summary

Properties

Form	Liquid
Storage instructions	Store at +4°C.
Storage buffer	Constituent: 100% Tween

Applications

Our [Abpromise guarantee](#) covers the use of **ab128987** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
IHC-P		Use at an assay dependent concentration.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors