

Product datasheet

Anti-WASH complex subunit 4 antibody ab111109

★★★★★ 1 Abreviews 1 References

Overview

Product name	Anti-WASH complex subunit 4 antibody
Description	Rabbit polyclonal to WASH complex subunit 4
Host species	Rabbit
Specificity	Previously labeled as WASH complex subunit 7
Species reactivity	<p>Reacts with: Human</p> <p>Predicted to work with: Mouse, Rat, a wide range of other species </p>
Immunogen	<p>Synthetic peptide corresponding to Human WASH complex subunit 4 (internal sequence). Database link: Q2M389</p>
General notes	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
Storage buffer	Constituent: Whole serum
Purity	Whole antiserum
Clonality	Polyclonal
Isotype	IgG

Target

Function	Component of the WASH complex, a complex present at the surface of endosomes that recruits and activates the Arp2/3 complex to induce actin polymerization. The WASH complex plays a key
-----------------	--

role in the fission of tubules that serve as transport intermediates during endosome sorting.

Sequence similarities

Belongs to the WASHS7 family.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors