

## Product datasheet

# Western Blot Stripping Buffer ab270550

1 Image

### Overview

<b>Product name</b>	Western Blot Stripping Buffer
<b>Tested applications</b>	<b>Suitable for:</b> WB
<b>General notes</b>	<p>Western Blot Stripping Buffer (ab270550) is a robust, ready-to-use solution allowing primary and secondary antibodies to be quickly removed from nitrocellulose and PVDF membranes and, as a result, Western blots to be reprobed with different antibodies. This permits more time efficient experiments that use less sample by reusing membrane without having to rerun gels and blots. Using a PVDF membrane will minimize loss of sample protein when stripping.</p> <p>This product is manufactured by Expedeon, an Abcam company. Expedeon product code ECLB0500 was previously called Western Blot Stripping Buffer - 500ml and is the same as the 500 mL size of this product. Expedeon product code ECLB5000 was previously called Western Blot Stripping Buffer - 5L and is the same as the 5000 mL size.</p> <p>The 500mL size is provided as 2 x 250 mL bottles.</p>

### Properties


<b>Form</b>	Liquid
<b>Storage instructions</b>	Shipped at Room Temperature. Store at Room Temperature. Store at +4°C.
<b>Storage buffer</b>	Constituent: 1.5% Hydrochloric acid

### Applications

**The Abpromise guarantee** Our [Abpromise guarantee](#) covers the use of ab270550 in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		Use at an assay dependent concentration.

### Images



Stripping and reprobing tests.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors