

Product datasheet

Anti-White Spot Syndrome Virus VP664 antibody ab26811

Overview

| | |
|----------------------------|--|
| Product name | Anti-White Spot Syndrome Virus VP664 antibody |
| Description | Rabbit polyclonal to White Spot Syndrome Virus VP664 |
| Host species | Rabbit |
| Tested applications | Suitable for: ELISA |
| Species reactivity | Reacts with: Other species |
| Immunogen | Synthetic peptide located on the putative domain of White spot syndrome virus VP664. |
| Positive control | synthetic WSSV VP664 peptideantigen. |

Properties

| | |
|-----------------------------|--|
| Form | Liquid |
| Storage instructions | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. |
| Purity | Immunogen affinity purified |
| Clonality | Polyclonal |
| Isotype | IgG |

Applications

Our [Abpromise guarantee](#) covers the use of **ab26811** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| Application | Abreviews | Notes |
|-------------|-----------|---------|
| ELISA | | 1/1000. |

Target

Relevance White spot syndrome virus [WSSV] is one of the most devastating shrimp pathogens , and it has

caused serious damage to the shrimp culture industry. This virus can infect several crustacean species, including shrimp, crab, and crayfish. VP664, is one of the viral structural proteins. The existence of the corresponding 18-kb transcript was confirmed by sequencing analysis of reverse transcription-PCR products, which also showed that vp664 was intron-less.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors