


Product datasheet

Anti-ZFY antibody ab155318

1 Image

Overview

Product name	Anti-ZFY antibody
Description	Rabbit polyclonal to ZFY
Host species	Rabbit
Specificity	ab155318 may cross-react with ZFX based on homology search showing these proteins (ZFX and ZFY) are 96% identical.
Tested applications	Suitable for: WB
Species reactivity	Reacts with: Human Predicted to work with: Cow 
Immunogen	Recombinant fragment corresponding to a region within amino acids 1 - 230 of Human ZFY (Uniprot ID P08048)
Positive control	NT2D1, U87-MG and SK-N-SH cell lysates

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
Storage buffer	pH: 7.00 Preservative: 0.01% Thimerosal (merthiolate) Constituents: 78% PBS, 20% Glycerol, 1% BSA
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab155318** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/1000 - 1/10000. Predicted molecular weight: 91 kDa.

Target

Function	Probable transcriptional activator.
Sequence similarities	Belongs to the krueppel C2H2-type zinc-finger protein family. ZFX/ZFY subfamily. Contains 13 C2H2-type zinc fingers.
Cellular localization	Nucleus.

Images



All lanes : Anti-ZFY antibody (ab155318) at 1/5000 dilution

Lane 1 : NT2D1 cell lysate

Lane 2 : U87-MG cell lysate

Lane 3 : SK-N-SH cell lysate

Lysates/proteins at 30 µg per lane.

Predicted band size: 91 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors