

Product datasheet

Anti-ZNF624 antibody - N-terminal ab228914

2 Images

Overview

Product name	Anti-ZNF624 antibody - N-terminal
Description	Rabbit polyclonal to ZNF624 - N-terminal
Host species	Rabbit
Tested applications	Suitable for: WB, ICC/IF
Species reactivity	Reacts with: Human
Immunogen	Recombinant fragment within Human ZNF624 (N terminal). The exact sequence is proprietary. Database link: Q9P2J8
Positive control	WB: SK-N-SH whole cell and nuclear extracts. ICC/IF: SK-N-SH cells.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
Storage buffer	pH: 7.00 Preservative: 0.01% Thimerosal (merthiolate) Constituents: PBS, 20% Glycerol
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab228914** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


Application	Abreviews	Notes
WB		1/500 - 1/3000. Predicted molecular weight: 100 kDa.
ICC/IF		1/100 - 1/1000.

Target

Relevance ZNF624 belongs to the krueppel C2H2-type zinc-finger protein family. It contains twenty one C2H2-type zinc fingers and may be involved in transcriptional regulation.

Cellular localization Nucleus

Images



All lanes : Anti-ZNF624 antibody - N-terminal (ab228914) at 1/500 dilution


Lane 1 : SK-N-SH (human neuroblastoma cell line) whole cell extract

Lane 2 : SK-N-SH (human neuroblastoma cell line) nuclear extract

Lysates/proteins at 30 µg per lane.

Predicted band size: 100 kDa

7.5% SDS-PAGE gel.



4% paraformaldehyde-fixed SK-N-SH (human neuroblastoma cell line) cells stained for ZNF624 (green) using ab228914 at 1/500 dilution in ICC/IF.

Blue: Hoechst 33342 staining.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors