



Product datasheet

Anti-ZNF674 antibody ab221463

1 Image

Overview

Product name	Anti-ZNF674 antibody
Description	Rabbit polyclonal to ZNF674
Host species	Rabbit
Tested applications	Suitable for: ICC/IF
Species reactivity	Reacts with: Human
Immunogen	Recombinant fragment corresponding to Human ZNF674 aa 89-121. ENSEMBL: ENSG00000251192 Sequence: VDEQIDHYKESQDKLPWQAAFIGKETLKDESGQ Database link: Q2M3X9  Run BLAST with  Run BLAST with
Positive control	U-2 OS cells.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
Storage buffer	pH: 7.2 Preservative: 0.02% Sodium azide Constituents: 40% Glycerol, 59% PBS
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab221463** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
ICC/IF		Use a concentration of 1 - 4 µg/ml. Fixation/Permeabilization: PFA/Triton X-100.

Target

Function	May be involved in transcriptional regulation.
Tissue specificity	Expressed in testis.
Sequence similarities	Belongs to the krueppel C2H2-type zinc-finger protein family. Contains 11 C2H2-type zinc fingers. Contains 1 KRAB domain.
Developmental stage	Expressed in fetal brain.
Cellular localization	Nucleus.

Images



Immunofluorescent analysis of PFA-fixed, Triton X-100 permeabilized U-2 OS cells labeling ZNF674 with ab221463 at 4 µg/ml. Antibody staining is shown in green.

Immunocytochemistry/ Immunofluorescence - Anti-ZNF674 antibody (ab221463)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors