

Product datasheet

Anti-ZNF692 antibody ab204595

2 References 2 Images

Overview

Product name	Anti-ZNF692 antibody
Description	Rabbit polyclonal to ZNF692
Host species	Rabbit
Tested applications	Suitable for: WB, ICC/IF
Species reactivity	Reacts with: Human
Immunogen	Recombinant fragment corresponding to Human ZNF692 aa 109-185. Sequence: GLVWECSAGHTFSWGPSLSPTPSEAPKPASLPHTTR RSWCSEATSGQELA DLESEHDERTQEARLPRRVGPPPETFP Database link: Q9BU19 Run BLAST with Run BLAST with
Positive control	RT-4, U-251 MG, Human Plasma and Human Liver lysats; HeLa cells.

Properties

Form	Liquid
Storage instructions	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
Storage buffer	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol
Purity	Immunogen affinity purified
Clonality	Polyclonal
Isotype	IgG

Applications

Our [Abpromise guarantee](#) covers the use of **ab204595** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

Application	Abreviews	Notes
WB		1/100 - 1/250. Predicted molecular weight: 57 kDa.
ICC/IF		Use a concentration of 1 - 4 µg/ml. Use PFA/Triton X-100 Fixation/Permeabilization.

Target


Function	May be involved in transcriptional regulation.
Sequence similarities	Belongs to the krueppel C2H2-type zinc-finger protein family. Contains 5 C2H2-type zinc fingers.
Cellular localization	Nucleus.

Images



Immunofluorescent analysis of HeLa cells (PFA-fixed/Triton X-100 permeabilized) labeling ZNF692 with ab204595 at 4µg/ml (green).

Immunocytochemistry/ Immunofluorescence - Anti-ZNF692 antibody (ab204595)



Western blot - Anti-ZNF692 antibody (ab204595)

All lanes : Anti-ZNF692 antibody (ab204595) at 1/100 dilution

Lane 1 : RT-4 cell lysate

Lane 2 : U-251 MG lysate

Lane 3 : Human plasma lysate

Lane 4 : Human liver lysate

Predicted band size: 57 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.com/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors